

OMEGA TYPE ES-206 Experimental Set Up is a combined setup for the three experiments covered under experimental Set Ups OMEGA TYPE ES-203, ES-204 and ES-205, for the study of experiments using Post Office Box. The set up is absolutely self contained and requires no other apparatus. Practical experience on this set up carries great educative value for Science and Engineering Students.



OBJECT

01. To determine the specific resistance of the material of a wire using Post Office Box.
02. To determine the internal resistance of a primary cell (Leclanche Cell) by Mance's method using Post Office Box.
03. To determine the resistance of a Galvanometer by Kelvin's method using Post Office Box.

FEATURES

The Set up consists of the following :

01. Post Office Box (Dial Type) OMEGA TYPE POB-180.
 02. Galvanometer 50-0-50 on bakelite Stand OMEGA TYPE MO-65.
 03. Leclanche Cell.
 04. Decade Resistance Box, One dial in step of 100 ohms, total 1000 ohms OMEGA TYPE DRBC-115R.
 05. Resistance wire of two different gauges each of $\frac{1}{2}$ metre length.
 06. Micrometer Screw Gauge.
- * Weight : 4 Kg. (Approx.)
 - * Adequate no. of connecting wires, 100cm long.
 - * Strongly supported by detailed Operating Instructions, giving details of Object, Theory, Design procedures, Report Suggestions and Book References.

We are committed to the continuous development of our products, and therefore reserve the right to amend specifications without prior notice.

OMEGA ELECTRONICS

India Sales

B-28, Fateh Singh Scheme , Opp. Rajputana
Palace Sheraton, Jaipur-302006, INDIA
Ph: +141-2375647 / 2379223
Fax: +91-141-2204481

Email: omega@sancharnet.in;
info@omegaelectronics.net
Web: www.omegaelectronics.net

Mfg./ Export/ Training

28 E & F, Malviya Industrial Area,
Jaipur-302017, INDIA
Ph: +91-141-2751136 / Fax: +91-141-2751559
Mobile: +91-9414072336