



OMEGA TYPE ES-304 Experimental Set-Up has been designed specifically for verification of Hooke's Law.

The set-up is complete in all respects and requires no other apparatus. Practical experience on this set up carries great educative value for Science and Engineering Students.

OBJECT

To verify the Hooke's Law.

FEATURES

The Set-up consists of the followings :

01 Hook's Law Apparatus : Consisting of 12cm long mirror scale, spring, slotted weights with hanger & pointer, with complete stand.

02 Strongly supported by detailed Operating Instructions, giving details of Object, Theory, Design procedures, Report Suggestions and Book References.

OTHER APPARATUS REQUIRED

01 DIGITAL STOP CLOCK : With START/STOP operation by means of toggle switch & RESET by a push button switch. It has a range of 999.9 seconds with resolution of 0.1 seconds and accuracy of $\pm 0.01\%$ (Quartz controlled). Display is through 4 no's of 12.5mm bright Seven Segment Displays and working voltage of the unit is $230V \pm 10\%$ 50Hz.

We are committed to the continuous development of our products, and therefore reserve the right to amend specifications without prior notice.

OMEGA ELECTRONICS