

**DETERMINATION OF E.C.E. OF COPPER
AND REDUCTION FACTOR OF A
HELMHOLTZ GALVANOMETER
OMEGA TYPE ES-360**



OMEGA TYPE ES-360 Experimental Set-up has been designed specifically to determine E.C.E. of copper and reduction factor of a Helmholtz galvanometer. The set-up consists of Copper Voltmeter, Helmholtz Galvanometer, Battery eliminator, Reversing key etc.

The set-up is complete in all respect and requires no other apparatus. Practical experience on this set-up carries great educative value for Science and Engineering Students.

OBJECT

01 To determine Electro-Chemical Equivalent of copper and reduction factor of a Helmholtz Galvanometer.

FEATURES

The Complete Experimental Set-up consists of the following :

01 Copper Voltmeter.

02 Helmholtz Galvanometer: Having two coils exactly parallel to each other so connected that current passes around them in the same direction. Two coils increase the uniformity of the field. Complete with Pye type Bakelite compass box and with windings of 2,50 and 500 turns.

03 Battery Eliminator : 0-5V D.C. at 3A, continuously variable regulated and short circuit protected OMEGA TYPE BE-5/3.

04 Reversing Key.

05 Adequate no. of connecting wires,

06 Strongly supported by detailed Operating Instructions, giving details of Object, Theory, Design procedures, Report Suggestions and Book References.

OTHER APPARATUS REQUIRED :

01 Physical Balance with Weight Box

We are committed to the continuous development of our products, and therefore reserve the right to amend specifications without prior notice.

OMEGA ELECTRONICS

Works:
28E & F, Malviya Industrial Area,
Jaipur-302 017 (INDIA)
Phone: 0141-2751559

E-mail : info@omegaelectronics.net
: omegajaipur62@gmail.com

Marketing Division:
B-28, Fateh Singh Scheme, Opp. Rajputana
Palace Sheraton, Jaipur-302006 (INDIA)
Phone : 091-141-2375647, 2379223

www.omegaelectronics.net